| | |
|--|---|
| BioGenex | |
| <div>48810 Kato Road, Suite 100E & 200E, Fremont,CA 94538</div> <div> </div> | <div>EmergoEurope, Prinsessegracht 20 2514 AP The Hague,, TheNetherlands</div> <div> </div> |
| <div>Tel : +1(800)421-4149, Fax: +1(510)824-1490, support@biogenex.com</div> <div> </div> | |

SUBSTRATE PACKS

| Doc. No. 932-HK092E-4 Rev. No.R |
|--|
| Release Date: 10-Aug-2020 |
| ENGLISH |
| Intended Use |
| Substrate packs are intended for in vitro diagnostic use as chromogens for detection of specific antigen-antibody reactions in immunohistochemistry (IHC) and in situ hybridization (ISH). |
| Summary and Explanation |
| BioGenex provides several substrate pack options both for manual and automated systems for in situ detection using alkaline phosphatase and peroxidase enzymes. The kits are designed to reduce substrate incubation times and minimize exposure to chemical hazards. These substrate packs allow microscopic visualization of different cellular components depending on the marker used either by in situ hybridization (ISH) or immunohistochemistry (IHC). |
| Principles of the Procedure |
| The demonstration of targets of interest in tissues and cells by immunostaining is a two-step process involving first, the binding of a primary antibody to the antigen of interest, and second, the detection and visualization of bound antibody by one of a variety of enzyme based chromogenic systems. The type of chromogenic substrate depends upon the type of enzyme used, e.g. AEC and DAB can be used with Horseradish Peroxidase; Fast Red can be used with Alkaline Phosphatase; Levamisole can be used with alkaline phosphatase substrate in the place of peroxide block to block endogenous phosphatase staining. In situ hybridization involves binding of specific nucleic acid probe that is labeled with an hapten. The label is subsequently detected using a primary antibody to the label followed by detection and visualization as above. |
| DAB (diaminobenzidine) substrate offers the greatest sensitivity of all the Horseradish Peroxidase colorimetric chromogens. The insoluble, permanent brown precipitate has a high-contrast in photographs. In addition, the sensitivity can be enhanced by carrying out the reaction in the presence of nickel or cobalt chloride and/or by examining slides by reflection interference microscopy (10-100x sensitivity). AEC (aminoethyl carbazole) is a colorimetric substrate for Horseradish Peroxidase. The bright brick-red reaction product is insoluble in water, but soluble in alcohol and xylene. The AEC substrate is suitable for immunohistochemistry (IHC), in situ hybridization (ISH), and membrane blotting applications. For IHC and ISH, the AEC substrate is compatible with aqueous mounting media. |
| Fast Red (4-chloro-2-methyl-benzenediazonium salt) is a substrate for Alkaline Phosphatase and offers high sensitivity for light microscopic observations. The bright red dye precipitate produces maximal contrast with blue counterstains and reproduces well by color photomicrography. The reaction product is insoluble in water, but soluble in organic solvents like alcohol and xylene.The dye precipitate is inherently fluorescent and sections can be examined by fluorescence microscopy, which results in a 5x increase in signal-to-noise ratio. Fluorescence emission of Fast Red precipitate is 540-590nm with a peak at 562nm. Excitation max is observed at 550nm. Visualization by fluorescence microscopy can be performed with the rhodamine, fluorescein or AMCA filter systems. The fluorescence is very stable and is not subject to photobleaching. |
| New Fuchsin |
| New Fuchsin, also known as Magenta III, produces an intense red/fuschia precipitate in the presence of alkaline phoshatase. It is not soluble in xylenes and alcohols, therefore it may be permanently mounted. |

| Substrate packs and Reagents available | | |
|--|--------------------------------|-------------|
| Catalog No. | Substrate Packs | Format |
| HK092-5KE | AEC SUBSTRATE PACK KIT | 250 slides |
| HK129-5KE | AEC SUBSTRATE PACK KIT, CONC | 2500 slides |
| HK124-7K, 9K | LIQUID DAB CHROMOGEN | 4 ml, 10 ml |
| HK130-5KE | LIQUID DAB SUBSTRATE PACK KIT | 250 slides |
| HK542-XAKE | 2- COMPONENT DAB PACK | 1000 slides |
| HK180-5K | NAPHTHOL PHOSPHATE SUBSTRATE | 5 ml |
| HK180-7K | NAPHTHOL PHOSPHATE SUBSTRATE | 15 ml |
| HK182-5KE | FAST RED SUBSTRATE PACK KIT | 500 slides |
| HK183-5KE | NEW FUCHSIN SUBSTRATE PACK KIT | 400 slides |
| HK139-06KE | ONE STEP AEC SOLUTION | 6ml |
| HK139-50KE | ONE STEP AEC SOLUTION | 50 ml |

Reagents Required but Not Supplied

All reagents required for IHC and ISH are not provided. See antibody and detection kit datasheet for complete set of reagents required for immunohistochemistry (IHC) or in situ hybridization (ISH) procedures.
Storage and Handling
Store all reagents at 2-8°C. Do not use after expiration dates as indicated on the reagent labels.

Precautions

HK092-5KE, HK129-5KE, HK542-XAKE:  Rep 2;R61.  Xn;R20/21.  Xi;R36 = May cause harm to the unborn child.Harmful by inhalation and in contact with skin.Irritating to eyes. S25 S26 S38 S45 S53 S36/37/39 S60 P11 = Avoid contact with eyes.In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.In case of insufficient ventilation, wear suitable respiratory equipment.In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).Avoid exposure - obtain special instructions before use.Wear suitable protective clothing, gloves and eye/face protection.This material and its container must be disposed of as hazardous waste.

HK183-5KE:  T;R23-25.  C;R35.  Xi;R37 = Toxic by inhalation.Toxic in contact with skin.Causes severe burns.Irritating to respiratory system. S26 S38 S45 S24/25 S36/37/39 S60 = In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.In case of insufficient ventilation, wear suitable respiratory equipment.In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).Avoid exposure - obtain special instructions before use.Wear suitable protective clothing, gloves and eye/face protection.This material and its container must be disposed of as hazardous waste.

Substrate packs and Reagents available

case of contact with eyes, rinse immediately with plenty of water and seek medical advice.In case of insufficient ventilation, wear suitable respiratory equipment.In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).Avoid contact with skin and eyes.Wear suitable protective clothing, gloves and eye/face protection.This material and its container must be disposed of as hazardous waste.

Staining Protocol

Refer to detection system manual for relevant staining protocols.

Preparation of Substrates and Chromogens

For Horseradish Peroxidase Kits:

AEC Kits: Prepare substrate solution by adding 1 drop of chromogen (AEC) to 1 vial (2.5 ml) of substrate and mix well. Each 2.5 ml of working solution is sufficient for up to 20 slides. This solution remains stable at room temperature (20-26°C) up to 5 hours.

DAB Kits: (HK542-XAKE) Add 2 drops (~80ul) of Liquid DAB Chromogen to 1 ml ready-to-use Substrate Buffer. This solution remains stable at room temperature (20-26°C) up to 6 hours.

For Alkaline Phosphatase Kits: Add 1 tablet of Fast Red chromogen to 1 vial (5 ml) of substrate (naphthol phosphate). Shake well until dissolved. This solution remains stable at room temperature (20-26°C) up to 8 hours. Optional: levamisole (not included) may be added to the Fast Red chromogen to reduce background staining. It serves as an inhibitor of endogenous alkaline phosphatase in most tissues.

For New Fuchsin Kits: Combine 50 µl (1drop) of New Fuchsin Chromogen Solution and 50µl (1 drop) of New Fuchsin Activator Solution in the bottom of the mixing vial. This is a critical step for successful color development. Mix reagents thoroughly by repeated gentle pipetting. Empty contents of one vial of Tris Buffer into Chromogen-Activator mixture prepared above and mix well. Add 400 µl (8 drops) of Substrate Solution to the solution prepared above, and mix well. NOTE: If endogenous alkaline phosphatase is suspected, add levamisole (HK113-5K) at a concentration of 0.6 mg/ml to the Substrate at this time to inhibit the alkaline phosphatase activity. Prepare fresh substrate immediately prior to use.

Procedure for application of substrate solution:

Refer to the antibody and detection kit datasheets for complete staining protocols and washing steps. User should follow BioGenex recommendations and validate any other conditions.

| Substrate | BioGenex Recommendations |
|-------------|-----------------------------------|
| DAB | 5-10 minutes at room temperature |
| AEC | 5-10 minutes at room temperature |
| Fast Red | 15-20 minutes at room temperature |
| New Fuchsin | 15-40 minutes at room temperature |

One may choose to monitor the positive control slide (the desired antibody versus the known positive tissue) under the microscope until acceptable color intensity has been reached. Rinse well with deionized water or rinse buffer.

Quality Control

Refer to the appropriate detection system package inserts for guidance on general quality control procedures.

Trouble shooting

Refer to the troubleshooting section in the package inserts of BioGenex Super Sensitive Detection Systems (or other equivalent detection systems) for remedial actions on detection system related issues, or contact BioGenex Technical Service Department at (925) 275-0550 to report unusual staining.

Expected Results

Staining using the IHC and ISH systems should result in deposition of colored chromogen pigment at the site of specific interaction with minimal to no non-specific background.

For AEC Kits: AEC (3-amino-9-ethylcarbazole) forms a brick red end product.

For DAB Kits: DAB (3,3'-diaminobenzidine) forms a brownish end product.

For Alkaline Phosphatase Kits: Naphthol Phosphate/Fast Red forms an intense red color.

For New Fuchsin Kits: New Fuchsin produces an intense red/fuschia precipitate.

Limitations of the procedure

Super Sensitive™ Detection kits demonstrate antigens that survive tissue fixation, embedding and sectioning. Correct treatment of tissues prior to fixation and embedding, while less critical for BioGenex Super Sensitive™ Reagents, is still important for obtaining optimal results. Inconsistent results may be due to variation in fixation and embedding methods employed by different laboratories, as well as from inherent variations in tissue. The results from immunostaining must be correlated with other laboratory findings and the relevant controls. An internal tissue processing control (e.g. vimentin) may be used to reveal errors in tissue processing. Use of BioGenex Antigen Retrieval pretreatment technique may permit recovery of antigenicity in formalin-fixed tissue. Please call BioGenex for more information on these products and their use in the standardization of immunostaining results.

The clinical interpretation of any positive staining or its absence should be evaluated within the context of clinical presentation, morphology and other histopathological criteria. The clinical interpretation of any positive staining or its absence should be complemented by morphological studies using proper positive and negative internal and external controls as well as other diagnostic tests. It is the responsibility of a qualified pathologist who is familiar with the proper use of IHC antibodies, reagents and methods to interpret the stained preparation. Staining is to be performed in a certified licensed laboratory under the supervision of a pathologist who is responsible for reviewing the stained slides and assuring the adequacy of positive and negative controls.

Tissues from persons infected with hepatitis B virus and containing hepatitis B surface antigen (HBsAg) may exhibit nonspecific staining with horseradish peroxidase. (Omata, et al. 1980)

Normal/ non-immune sera from the same animal source as secondary antisera used in blocking steps may cause false-negative or false-positive results due to autoantibodies or natural antibodies.

False-positive results may be seen due to non-immunological binding of proteins or substrate reaction products. They may also be caused by pseudoperoxidase activity (erythrocytes), endogenous peroxidase activity (cytochrome C), or endogenous biotin (e.g., liver, breast, brain, kidney) depending on the type of immunostain used. (Nadji & Morales, 1983)

Performance Characteristics

BioGenex has conducted studies to evaluate the performance of all its substrate packs using several BioGenex IHC and ISH assays. BioGenex substrate packs have shown reproducible and consistent results when used within a single run, between runs, between lots and wherever applicable between manual and automated runs. The products have been determined to be stable for the periods specified on the labels either by standard real time or accelerated testing methods. BioGenex ensures product quality through 100% quality control for all products released and through surveillance programs.

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